











Symbol	Title	Description	Reference
	Manufacturer	Indicates the medical device manufacturer	ISO 15223-1 5.1.1
	Use-by date	Indicates the date after which the medical device is not to be used	ISO 15223-1 5.1.4
	Batch code	Indicates the manufacturer's batch code so that the batch or lot can be identified	ISO 15223-1 5.1.5
	Catalog number	Indicates the manufacturer's catalog number so the device can be identified	ISO 15223-1 5.1.6
	Country of manufacture	To identify the country of manufacture of products	ISO 15223-1 5.1.11
	Sterilized using irradiation	Indicates a medical device that has been sterilized using irradiation	ISO 15223-1 5.2.4
	Do not use if package is damaged	Indicates that a medical device should not be used if the package has been damaged or opened	ISO 15223-1 5.2.8
	Do not re-use	Indicates a medical device that is intended for single use only	ISO 15223-1 5.4.2
	Caution	Indicates that caution is necessary when operating the device or control close to where the symbol is placed, or that the current situation needs operator awareness or operator action in order to avoid undesirable consequences	ISO 15223-1 5.4.4
	Unique device identifier	Indicates a carrier that contains unique device identifier information	ISO 15223-1 5.7.10